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TABLE OF CONTENTS

Part 1. Technical Passport

1. Safety rules.....	68
2. Intetnded purpose.....	73
3. Unit description.....	74
4. Procedure instruction.....	78
5. Package.....	86
6. Maintenance and battery replacement procedure.....	87
7. Troubleshooting.....	90
8. Technical information.....	93
9. Manufacturer`s guarantees.....	101
10. Manufacturer`s address.....	103

Part 2. Application Instruction

1. General information.....	106
-----------------------------	-----

2. Indications and contraindications.....	107
3. Treatment conditions.....	111
4. Electrostimulation intensity.....	113
5. Stimulation methods.....	115
6. Indication of modes.....	116
7. Recomendation on selection of treatment areas.....	120

Annex 1. Segmental dermatomeres

(segmental zones).....	185
------------------------	-----

Annex 2. Accessories used.....199

Warranty Maintenance Form.....	205
--------------------------------	-----

Acceptance certificate.....	208
-----------------------------	-----

Thank you for using the DENAS unit.

In order to make the usage of the unit effective and safe, please read carefully all section of this manual.

PART 1

Technical Passport

1. SAFETY RULES

Please read carefully all information marked with this sign. This information is important for ensuring safety and effective use of the unit.

The information contained in this Operating Manual is important for your safety and for correct application and maintenance of the unit.

The unit is safe in operation as it uses the low voltage power supply isolated from the working part of the unit (type B device with a type F working part).

The unit may not be used to treat patients with implanted electronic devices (such as the heart pacemaker) and to treat patients with idiosyncrasy to electric current.

The unit may not be used in the area of direct anterior projection of the heart.

During stimulation the patient should not be connected to any high frequency electric device. Simultaneous use of the unit and other electric devices may cause burns and possible damage of the unit.

Operation in the vicinity of short-wave or microwave equipment may cause the output parameters of the unit to become unstable.

The product contains fragile components. Please avoid shocks.


Avoid long-term exposure of the unit to direct sunrays at a high ambient temperature (over 25 °C). This may cause failure of the indicator on the unit or destroy its batteries.


The unit is not waterproof. Protect it from moisture.


All maintenance works on the unit must be performed by skilled technicians at the manufacturing factory.

Transportation conditions: temperature -50 to +50 °C, relative air humidity 30 to 93 %, atmospheric pressure 70 to 106 kPa (525 to 795 mm Hg).

Storage conditions: temperature -50 to +40 °C, relative air humidity 30 to 93 %, atmospheric pressure 70 to 106 kPa (525 to 795 mm Hg).

Operating conditions: temperature +10 to +35 °C, relative air humidity 30 to 93 %, atmospheric pressure 70 to 106 kPa (525 to 795 mm Hg).

Warning! *If the unit was stored at the ambient temperature lower than 10°C, before operating keep it under normal climatic conditions for at least two hours.*

Recycling: the unit packaging materials are not hazardous for the environment; they can be recycled.

Separate collection of electric and electronic equipment.

Warning! *The unit contains precious materials that may be used repeatedly after disposal with consideration of the environmental protection regulations. They must be submitted to specially designated collection and recycling facilities (please consult in the relevant services in your locality).*

Warning! *No modification of this equipment is allowed.*

2. INTENDED PURPOSE

The DENAS unit is used for electrostimulation of biologically active points and areas. The unit is equipped with integrated electrodes and a plug for connecting detachable electrodes.

The DENAS unit is designed for medical and preventive treatment facilities and individual application.


Compliance:

This medical device bears the CE marking in compliance with Directive 93/42/EEC concerning Medical Devices.

3. UNIT DESCRIPTION


EN

Picture 1. Arrangement of the unit. Front



***Warning!** This plug is only designed for connection of electrodes manufactured by the supplier. Never connect the unit to the computer as it may cause failure of the computer and expose the user to the risk of electrocution.

Picture 2. Arrangement of the unit. Back

* The unit contains a duplicate label with the serial number under the batteries (inside the battery compartment) of the unit. Make sure not to remove it and not to delete the serial number as it results in loss of warranty for the unit.

Picture 4. Designation of symbols on the unit display

4. PROCEDURE INSTRUCTION

4.1. Procedure conditions

No special conditions are required for the procedure. The room should be dry and warm. During the procedure the patient should be positioned comfortably, seated or lying. After the procedure it is recommended to remain relaxed for 10 – 15 minutes.

4.2. Switching on the unit

Press the  button; the unit should turn on and switch to the standby mode.

4.3. Mode selection

To select the desired mode, use the  (Mode) button.

Each time the  button is pressed, the modes are switched on in the following

sequence: Test → 10 Hz → 77 Hz → 140 Hz ( →  → ).

4.4. Stimulation power setting

In order to set the desired power, use the  and  buttons:

- each time the button is pressed, the power is decreased or increased by one unit;
- for continuous increase or decrease of the stimulation power hold down the button.
When the desired power level is reached, release the button.

The set stimulation power level is shown in the left top corner of the display as a number symbol and on the linear horizontal scale.

4.5. Timer

If you need to set a fixed stimulation time, use the Timer function intended for operation in the Therapy 1 and (77 Hz) and Therapy 2 (140 Hz) modes and is switched on by pressing the  (Timer) button.

When the Timer function is switched of, the maximum duration of the procedure is 20 minutes.

Each time you press the button, the stimulation time is set at 5, 10, 15 and 20 minutes and the Timer on indication is displayed on the unit.

Next pressure of the button (after a 20-minute interval) switches off the Timer function; the Timer On indication on the display disappears.

4.6. Volume setting

The unit can operate in a mute mode or you can select one of the two sound alarm volume settings. The volume level is set using a combination of the ,  and  buttons.

To select the sound alarm volume, press a combination of  and  buttons simultaneously to increase the volume or the  and .

The volume level is indicated on the display of the unit: no picture; ;).

4.7. Activation of the Forced Therapy mode.

In the Forced Therapy mode, the unit can operate in the therapy modes irrespective of the presence and the quality of the contact of the electrodes with the skin surface.

To activate the Forced Therapy function, press simultaneously the and buttons; the mode is deactivated by simultaneously pressing the and buttons.

Warning! *With the Forced Therapy Function activated, the unit is not switched off automatically when there is no contact between the electrodes and the skin, even if the keyboard buttons are not used. Do not leave the unit with the Forced Therapy function activated as it will cause discharge of the battery before the indicated time.*

4.8. Switching off the unit

Press the  button and hold it pressed until the unit switches off.

4.9. Preparation of the unit for the procedure.

Before using the unit, make sure that there are no impact marks on its case. The liquid crystal display of the unit is not damaged. Make sure that there are no traces of foreign materials or marks of damage on the electrode.

Warning! Do not use the unit if it is damaged.

Warning! If the unit was stored in a cold place or a place with high humidity, keep it at room temperature for at least 2 hours.

4.9.1. Open the consumer package and take out the unit.

4.9.2. If there are no batteries in the battery compartment, install them in accordance with paragraph 6.3 of this manual.

4.10. Procedure

4.10.1. Prepare the skin and the surface of the electrodes for the procedure.

4.10.2. Switch on the unit, using the  button.

4.10.3. If necessary, set the sound alarm volume as specified above and fixed duration of the procedure, using the Timer function.

4.10.4. If required, connect one of the detachable electrodes to the respective plug (picture 1) specified in Annex 2. In case of a small skin contact area or in case of a dry skin turn on the Forced Therapy function on the unit as per paragraph 4.7.

Carry out the procedure as specified below. Upon completion of the

procedure turn on the Forced Therapy function.

EN

Warning! Do not connect to the unit detachable electrodes not provided by the manufacturer.

4.10.5. Set the required stimulation mode, using the  button.

4.10.6. Before the procedure place the electrodes on the skin in the stimulation area, making sure that they are in good contact with the skin. The contact symbol  will be displayed.

4.10.7. Gradually increase the power until it reaches the desired level. As soon as the power level is above 1, the stimulation time will start.

Warning! It is not recommended to carry out the procedure at an intolerable power level!

If the Timer function deactivated, the procedure time count in the order of increase starts on the digital indication display. If the Timer function is activated, the time count on the digital indication display starts in the order of decrease from the set timer value.

- 4.10.8.** When the maximum procedure time (30 minutes) is reached or earlier (if the timer is activated), the alarm sounds and electrostimulation stops. To continue stimulation, either remove the unit from the skin for a short time or change the electrostimulation power or mode. The alarm also sounds when the electrode loses contact with the skin.

Warning! *Inadequate cleaning of the electrodes may lead to allergic reactions or cross-infection in case of use of the unit by several patients.*

Warning! *In case of idiosyncrasy, in case of aggravation of the patient's condition during the procedure or skin irritation after the procedure, stop the stimulation and consult a physician.*

5. PACKAGE

The package contains the following items:

Name	Number, pieces
DENAS unit	1
Operation manual	1
Type LR6/AA battery	2*
Consumer packaging	1

* *The number of batteries may be increased.*

6. MAINTENANCE AND BATTERY REPLACEMENT PROCEDURE

6.1. Daily maintenance should include the following operations:


6.1.1. External inspection of the unit. Make sure that there are no impact marks on the unit as this may lead to incorrect operation of the unit.

6.1.2. Clean the integrated electrode and the detachable electrodes (if they were applied) before and after use. The lack of adequate cleaning may cause allergic reactions or skin infection, especially in case of use of the unit by several patients. Use standard disinfectants (such as 70% alcohol solution) and soft pile-free tissues.

6.2. If the unit will not be used for a long period of time (more than 3 days), remove the batteries from the battery compartment. This will prevent their early discharge.

6.3. If the battery indicator on the LCD () starts flashing, replace the batteries. To do this:


- open the battery compartment as shown in picture 5;
- remove the used batteries;
- install new batteries*, making sure that the polarity is correct;
- close the battery compartment cover on the unit.

1. Remove the battery compartment

2. Install new batteries, making sure that the polarity is correct

3. Close the battery compartment cover

Picture 5. Unit battery replacement procedure

7. TROUBLESHOOTING

The table below describes the situations that may be interpreted as fault and can be remedied without external help. In case of other defects contact the manufacturer's representative (see Section 9 "Manufacturer's Warranty"), do not try to remedy them. Possible faults and remedies.

Problem	Possible cause	Remedy
The unit does not turn on when the  button is pressed.	No batteries	Install new batteries (see Maintenance and Battery Replacement Procedure)
The unit does not turn on when the  button is pressed.	Battery voltage less than 2.1 V	Replace batteries (see Maintenance and Battery Replacement Procedure)

Problem	Possible cause	Remedy
When the unit is turned on, alarm sounds and the unit switches off automatically	Battery voltage less than 2.2 V	Replace batteries (see Maintenance and Battery Replacement Procedure)
	Time interval between de-energizing and consequent energizing of the unit less than 3 sec.	Keep the unit de-energized for at least 3 sec. and then turn it on again
The unit switches off automatically or the battery symbol flashes and intermittent alarm sounds	Battery voltage less than 2.1 V	Replace batteries (see Battery Replacement Procedure)

Problem	Possible cause	Remedy
The unit does not switch to the stimulation mode when the integrated electrodes are in contact with the skin surface; no skin contact symbol on the display	Contact area between the integrated electrodes and the skin not sufficient	Apply the integrated electrodes of the unit to the skin unit, making sure the contact is sufficient. If necessary, moisten the contact area or turn on the Forced Therapy function, using a combination of the Power  and  buttons.
No sound alarm	The sound is off	Turn on the sound, using the  and Power  buttons
Quick consumption of batteries	Poor quality batteries	Use good quality batteries (alkaline batteries are recommended) or rechargeable batteries of the corresponding type and size with the voltage not exceeding 1.5 V
	Forced Therapy function activated	Deactivate the Forced Therapy function, using a combination of the Power  and  buttons.

8. TECHNICAL INFORMATION

8.1. Maximum consumption current of the unit at the supply voltage of $3,0 \pm 0,1$ V;
— in the stimulation mode – max. 140 mA;
— in a de-energized condition – max. 10 μ A;

8.2. Weight without batteries – max 300 g.

8.3. Dimensions – max. 138×52×37 mm.


8.4. The unit turns off automatically if not used for more than 3 min. (except in the Forced Therapy mode).

8.5. Operating parameters:

8.5.1. Stimulation impulse frequency: Test mode — 10 ± 1 Hz; Therapy 1 mode — 77 ± 7 Hz; Therapy 2 mode 140 ± 0.5 Hz

8.5.2. Parameters of individual stimulation impulses:

8.5.2.1. Individual impulse waveform with applied load of $(20 \pm 10\%)$ kOhm:



8.5.2.2. At minimum power, load $(20 \pm 10\%)$ kOhm:

- amplitude in the positive part – max 30 V;
- amplitude in the negative part – max 10 V.



8.5.2.3. At maximum power, load ($20 \pm 10\%$) kOhm:

- amplitude in the positive part – (30 ± 10) V;
- amplitude in the negative part – (300 ± 70) V.



8.6. Relation between the impulse waveform and the load resistance at the medium and maximum power level:


Load resistance	Power level: Half		Power level: Full	
	Waveform	V_{p-p}	Waveform	V_{p-p}
Free output		$\approx 280 \text{ V}$ $\pm 20 \%$		$\approx 330 \text{ V}$ $\pm 20 \%$


200 Ohm

 $\approx 15 \text{ V} \pm 20 \%$

 $\approx 15 \text{ V} \pm 20 \%$

500 Ohm

 $\approx 20 \pm 20\%$
 $E_{\text{pulse}} \approx 48 \text{ uJ}$

 $\approx 35 \text{ V} \pm 20\%$
 $E_{\text{pulse}} \approx 170 \text{ uJ}$

8.7. Effective current – I_{eff} for load resistance 500 Ohm:

Mode	Power level: Half	Power level: Full
1Th (77 Hz)	$I_{\text{eff}}=2,7 \text{ mA}$	$I_{\text{eff}}=5,1 \text{ mA}$
2Th (140 Hz)	$I_{\text{eff}}=3.7 \text{ mA}$	$I_{\text{eff}}=6.9 \text{ mA}$

8.8. Electromagnetic radiation

Test	Conformity	Conditions of use
HF radiation CISPR 11	Class B	The DENAS unit can be used in all kinds of facilities, including application at home

8.9. HF radiation resistance

Test	IEC 60601-1-2, testing conditions	Acceptable level
IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms
IEC 61000-4-3	3 V/m, 80 MHz to 2,5 GHz	3 V/m

8.10. Electromagnetic field resistance

Test	Test level	Conformity level	Conditions of use
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 4 kB contact ± 8 kB air	The floor should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic materials, relative air humidity should be at least 40 %

Magnetic fields IEC 61000-4-8	3 A/m	3 A/m	The magnetic field parameters should not exceed the limits for commercial buildings and medical facilities
----------------------------------	-------	-------	--

8.11. Recommendations concerning determination of the required distance between the DENAS stimulation unit and radio frequency radiating equipment.

Declared maximum output power of the transmitter P (W)	Radiation frequency and the d distance determination formula (m)		
	150 kHz — 80 MHz $d = 1,2\sqrt{P}$	150 kHz — 800 MHz $d = 1,2\sqrt{P}$	800 MHz — 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

9. MANUFACTURER'S GUARANTEES

9.1. Lifetime of the unit — 5 years.

If the operating instructions are followed, the lifetime may be much longer than the officially established one.

9.2. Warranty operation term — 24 months of the date of sale.

9.3. The Seller (Manufacturer) or the organization acting on behalf of the Seller (Manufacturer) under the contract with him shall not be liable for defects if they have emerged after the transfer of the unit to the user due to:

- 1) failure of the user to observe the transportation, storage, maintenance and operation rules provided for by this Manual;
- 2) mechanical damage;
- 3) third parties' actions;
- 4) circumstances of force majeure.

9.4. The warranty does not cover the units with a destroyed warranty label in the battery compartment.

9.5. In case of failure of the unit or its defect during the effective term of the warranty and if the owner discovers that certain components are missing, the owner of the unit shall send to the manufacturer or his representative a repair (replacement) request stating his full name, address and telephone number and containing a short description of the defect and the conditions in which and the date when it was discovered.

10. MANUFACTURER'S ADDRESS

RC ART LLC

Academician Postovsky St. 15, Ekaterinburg 620146, Russia

Tel.: +7 (343) 267-23-30

<http://www.denascorp.ru>

e-mail: corp@denascorp.ru

Official Representative in the European Union:

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Application Instruction

1. GENERAL INFORMATION

Stimulation of reflexogenous zones and points with impulse currents for preventive treatment, therapy and recovery of bodily functions is an independent method of treatment and can be applied both as part of comprehensive treatment and as individual therapy. Numerous studies show that electroneurostimulation (DENS) is based on multi-level reflectory and neurochemical reactions that trigger a cascade of regulatory and adaptive mechanisms of the body, which leads to elimination of pain syndromes and blood circulation improvement, has an anti-inflammatory effect, activates production of biologically active substances and metabolism in the tissues and normalizes muscular and vascular tonus. Dynamic electroneurostimulation promotes general health improvement and fitness.

2. INDICATIONS AND CONTRAINDICATIONS

The unit may be applied for treatment of patients any age from newly born children to elderly people.

Indications:

- acute and chronic pain syndromes;
- injuries;
- treatment of functional disorders in case of diseases of the respiratory apparatus, ENT organs, the cardiovascular, musculoskeletal, nervous, endocrine and urinogenital systems, the eyes and the skin;
- rehabilitation after diseases, surgery and injuries;
- enhancement of adaptability of the body in case of exposure to adverse factors and in the conditions of hard physical and intellectual work, in case of physical and intellectual overstress, chronic fatigue syndrome, difficulties waking up in the morning and sleepiness during the day, in case of dyskoimesis in the evening

and insomnia, increased irritability and depression and for prevention of catarrhal conditions.

Treatment effect

When using the unit, one should remember that the most apparent therapeutic effect is achieved after a course of treatment. Its duration is 6 to 8 procedures for some diseases, 8 to 12 procedures for other diseases or in some cases – 14 to 20 procedures. The combined therapeutic effect of individual procedures also creates an after-effect of the treatment course that continues after its completion and the long-term results are in some cases even more favourable than the direct ones. But one should not forget that long-term continuous exposure to one physical factor may cause the body to adapt to it, which may considerably decrease the treatment effect.

Moreover, the stochastic nature of physiological processes leads to the inadequate effect of or the lack of favourable response to therapeutic physical factors in some patients.

Warning! Application of the unit does not rule out the use of other treatment methods indicated by the attending physician.

Warning! In case of lack of positive changes or aggravation of the patient's condition consult the attending physician.

Contraindications

Absolute:

- idiosyncrasy to electric current;
- an implanted heart pacesetter;

Relative:*

- status epilepticus;
- neoplasia of any etiology and localization;

* In cases listed above consent of the physician is needed before indicating the treatment.

- acute fever of unknown etiology;
- venous thrombosis;
- acute mental, alcoholic or narcotic excitement

Warning! *If you have any of the above listed contraindications, consult your physician before using this device.*

Warning! *In case of idiosyncrasy, if you feel worse or the pain increases during the procedure, stop the stimulation and consult a physician immediately.*

3. TREATMENT CONDITIONS

Warning! Sudden pain of any localization may be the first and frequently the only sign of a serious disease. That is why if pain attacks happen for the first time and repeat and the pain becomes more intense, consult the physician as soon as possible.

DENS procedures do not require any special conditions. The procedures can be performed without attendance of a physician or a nurse or by an operator who will carry out the treatment in the zones and points that the patient cannot reach without help. During the procedure, the patient may be seated in an armchair or lie in a comfortable position. It is recommended to take a 10-15 minute rest after the procedure.

These instructions describes the stimulation schemes used in case of the most common diseases. Additional information on possible applications of DENS can be found in a manual on DENS therapy*.

* Manual on dynamic electrostimulation using DiaDENS and DiaDENS-DT units. – Ekaterinburg, 2005.

Warning! After each procedure clean the electrodes with a standard disinfectant solution (such as 70% alcohol solution). The unit should be stored with dry electrodes.

4. ELECTROSTIMULATION INTENSITY

Electrostimulation intensity (power) is set individually depending on the patient's individual sensations. For convenience electrostimulation power is divided into three levels (power ranges).

The first, or *minimum level* – you do not have any subjective sensations or feel a slight vibration in the area under the electrodes. This level is used when intensive stimulation is not recommended, that is for patients with low intensity pain. In case of cardiovascular diseases (arterial hypertension, arterial hypotension and vegeto-vascular dystonia syndrome) the procedure should also start with the stimulation at the minimum power level.

The second, or *comfortable level* – the patient feels vibration, pleasant tingling or a slight “burning” but without any painful sensations.

This level is used for treatment of patients complaining of slight and medium intensity pain. This is the most frequently used power level.

The third, or *maximum level* – the patient has a painful tingling or burning sensation. Stimulation of such intensity may be accompanied by involuntary contraction of the muscles located close to the electrodes.

This level is used for the treatment of a strongly manifested pain syndrome.

Warning! *At different stages of treatment the electrostimulation power level can be increased and decreased as the patient's individual tolerance changes and as the pain decreases.*

Warning! *The power level should be adjusted depending on the patient's sensations. Do not exceed the pain threshold.*

5. STIMULATION METHODS

During the procedure, you can use three stimulation methods — *the stable, the labile and the combined labile and stable method*.

The *stable method*, (fixed position of the electrodes) is applied for treatment of small areas and areas with skin changes (post-surgery and post-burn scars, oedemas etc.).

When using the *labile treatment method*, move the built-in electrodes smoothly around the treated area so that they remain in constant contact with the body surface, at the speed of 0.5 to 3 cm/sec. You may make linear, spiral, circular or other movements, depending on the area and relief of the treated zone, and only on undamaged skin.

To apply the *combined labile and stable treatment method*, move the electrodes across the skin surface, keeping the built-in electrodes in a fixed position in some points for some time (for example, in maximum pain areas). The pressure degree depends on the patient's subjective sensations.

6. INDICATION OF MODES

6.1. Therapy mode

Stimulation in the Therapy mode is used in case of local pain syndrome, functional disorders and as an emergency care method. Duration of the procedure in the Therapy mode in the area of direct projection of pain and functional disorders is selected based on the following response:

- the complaint is fully eliminated;
- the patient feels considerably better;
- bright redness of the skin or a feeling of warmth in the area under the electrodes.

The following frequencies are available in the Therapy mode:

- a) Therapy 1 — *therapeutic frequency 77 Hz*, is used in the direct projection area of the complaint and in segmental zones. The effect is achieved after 15–30 minutes and continues up to 1–1.5 hours and more.

Indications: functional disorders and pain syndromes in case of internal diseases and during the post-surgery period;

b) Therapy 2, *therapeutic frequency 140 Hz*, is used in the area of direct projection of the complaint. The effect is achieved in the course of the first few minutes and continues from several minutes to an hour. To enhance the treatment effect, you may then continue stimulation using the device at a lower frequency.

Indications: acute pain related to diseases and affection of the musculoskeletal system at the acute stage and the peripheral nervous system damage.


6.2. Test mode (preventive treatment)

Stimulation in the Test mode is used for treatment of any complaints to enhance effect of the procedure. Treatment can be carried out in the areas symmetrical to the complaint area and in the general (universal) action areas, such as the hand area (he-gu zone) as shown in picture 6.

The procedure should be carried out at the comfortable power level (slight tingling).

If you do not feel anything, press the  button several times until you feel a slight tingling. If you have unpleasant (pain) sensations under the electrode, press the  button until you feel a slight tingling.

The unit may be moved to the next area after the mode time is over (sound alarm or indication on the display).

Picture 6. Stimulation on the hand (he-gu)

7. RECOMMENDATIONS ON SELECTION OF TREATMENT AREAS

Warning! *Only skin can be exposed to stimulation by means of the electrodes.*

Direct stimulation of the lesion focus – in the lumbar area in case of lumbodysnia or on the skin in the area around the affected joints in case of pain in the knee joints – is one of the key, simple and rather effective methods of treatment area selection. To enhance the treatment effect, stimulation of the lesion focus can be combined with stimulation of the respective segmental reflexogenous areas (Annex 1).

It is recommended not to treat more than 2 or 3 zones during one procedure.

Recommended duration of one procedure:

- for children aged 0 to 12 months — 5-10 minutes;
- for children aged 12 to 36 months — 10-15 minutes;

- for children aged 3 to 5 years — 15-20 minutes;
- for children aged 5 to 12 years — 20-25 minutes;
- for children older than 12 years of age and adults — up to 30 minutes.

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Picture 1. Segmental innervation of the human body (anterior surface)

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Picture 2. Segmental innervation of the human body (posterior surface)

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Some indications for treatment of segmental dermatomes

Dermatomes	Symptoms, states and diseases
C1-C2	headache, pain in the throat and nose, pain in the occiput, pain in the shoulder area, hemiplegia;
C3-D1	headache, pain in the ears, pain in the throat, pain and tension of the occiput muscles, pains in the shoulder, back and loins
D1-D2	hemicrania and other headaches, tension and pain in the spine, pain in the shoulder blade, pain in the knee, pain in joints
D2-D3	pain and tension of muscles in the back, loins, shoulder, occiput area, intercostal neuralgia
D3-D4	pain and tension in the occiput area, pain in the exterior surface of the shoulder, shoulder blade, thorax area, in the loins, abdomen area, lumbodinia, pain in the sacrum and back

D4-D5	headache, pain in the neck area, pain the shoulder blade, in the area of shoulder and back, pain in area of the heart, pains in the stomach area
D5-D6	headache, pain in the back and thorax, intercostal neuralgia, pain in the spine
D6-D7	pain in the back and neck, restraint of the spine movement, intercostal neuralgia, pain in the back and heart, pain in the stomach
D7-D8	pain in the loins, pain in the abdomen
D9-D10	pain in the loins, pain in the thorax and hypochondrium, pain in the heart area, pain in the stomach
D10-D11	pain in the abdomen, loins and back, intercostal neuralgia, cough with pain, pain and heaviness in the abdomen
D11-D12	pain in the back, pain in the stomach, pain in the intestines
D12-L1	pain in the abdomen, back and spine, pain in the pit of the stomach
L1-L2	pain and contractions in the lumbar area, tension of the spine and loins muscles, pains and spasms in the stomach

L2-L3	headache, pain and tension in the back and loins, sensation of tension in the spine muscles, pain in the hip, lumbodynia, intestinal colic (pain in the intestines), pain in the external genitalia, pain in the area of external genitalia
L3-L4	pain in the loins
L4-L5	pain in the loins and lateral part of the pelvis; pain in the external area of the knee joint, lumbodynia, intestine colic
L5-S1	lumbodynia, lumbago, ischias
S1-S2	pain in the abdomen, sacrum and hip joint, lumbago, intestine colic
S2-S3	pain in the sacrum, loins, spine, weakness in the knee joint, pain in the joints, labor pains
S3-S4	pain in the loins, ischias, pain in the spine, intestine colic
S4-S5	pain in the loins and back, sacral-coccygeal pains, pain in the area of lateral surface of the buttock, lumbodynia, ischias, cystodynia, pain in the urethra

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

*A set of detachable zonal electrodes **DENS-applicator**. Designed for non-invasive therapy of the pain area, lesion foci and reflexogenous areas (with undamaged skin).*

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Therapeutic detachable massage electrodes. Designed for stimulation in the area of the pain, lesion foci, reflexogenous areas and acupuncture points of the body.

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
*Detachable paraorbital electrode **DENAS-goggles**. Designed for preventive therapy and treatment of biologically active points around the eyes in the paraorbital area using the dynamic electroneurostimulation method.*

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Detachable therapeutic electrode. Designed for therapy of biologically active points using the dynamic electroneurostimulation method.

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